

Magnetic Mat









The magnetic mat is made from Silicone Material. It utilizes specially designed magnets, which retain their magnetism for an extended period without any noticeable decrease in magnetic strength. The high grade silicone material enables the mat to be cleaned and re-sterilised by conventional methods (steam).

Use

Holding surgical instruments safely and securely is always a top priority in the operating environment. The Magnetic Mat has been specifically designed and tested within operating theatres. The flexible, slip-resistant design allows the mat to be placed directly upon the patients body giving the surgeon secures and easy access to their instruments. Magnetic Mat can be placed singularly or overlapped onto trolleys and theatre tables to prevent instruments from rolling or moving. This protects instruments and helps to keep them in their correct position until they are needed, improving accessibility for the surgical assistant.

Product

The Magnetic Mat is made from surgical grade silicone rubber containing magnets specially designed for steam sterilisation. The positions of the magnets within the Magnetic Mat are indicated by moulded areas. This allows instruments to be placed directly onto a magnet to maximize holding power. The Magnetic Mat features raised nipples and a raised lip, these additional features prevent instruments moving or rolling off the mat. The moulded silicone structure of the Magnetic Mat follows the contours of any surface it is placed on whilst maintaining its holding capability.

The Magnetic Mat gives the user the assurance of holding instruments with the maximum amount of security and will help minimize damage to instruments used within the operating theatre.

Importer:

Health Care Inc.

Office: 208 & 310, Arjun Avenue, Opp. Samartheshwar Mahadev, Nr. Law Garden, Ellisbridge, Ahmedabad - 380 006 (GUJARAT) INDIA

Health Care Inc. Ph.: 91-79-26431051 Fax: 91-79-26568532 E-mail: healthcare_zu@yahoo.com