MEDICAL STOCKING

MANUFACTURERS SINCE 1993




HOW VARICOSE VEINS ARE FORMED


The blood is drained from the lower limbs by means of the deep and superficial venous systems. These two systems are connected to each other by means of perforating veins, the valves of which only allow the blood to pass from the superficial system into the deep system, and not the other way around. Varicose veins are the result of a pathological condition that's primarily characterized by the expansion of large superficial veins and perforating veins.

There can be various causes, including genetics, foot or posture defects, gender (varicose veins are more commonly found in women), pregnancy, obesity and occupations that involve standing or sitting for several hours at a time. Graduated compression elastic stockings represent an excellent means of preventing and treating varicose veins.

HEALTHY VEIN CIRCULATION



VARICOSE VEINS



increasing the pressure inside the vein itself: this dilation results in varicose veins PERFORATING VEIN

DEEP VENOUS SYSTEM (POPLITEAL VEIN)

SUPERFICIAL VENOUS SYSTEM (LARGE SAPHENOUS VEIN)

In normal veins the blood flows from the bottom upwards.

In varicose veins the blood flow slowly in the opposite direction or rather towards the feet.

PERFORATING VEINS

Distributor in India :



Office: 208 & 310, Arjun Avenue, Opp. Samartheshwar Mahadev, Nr. Law Garden, Ellisbridge, Ahmedabad · 380 006 (GUJARAT) INDIA
Ph.: 91-79-26431051 Fax: 91-79-26568532 M. +91 98250 42470 E-mail: healthcare_zu@yahoo.com

ELASTIC STOCKINGS



Art. H52 - AD - mm/Hg 18-20 (at the ankle) Antiembolism below knee stocking with inspection hole. Comp: 82% Polyamide, 18% Elastane



Art. H51 - AG - mm/Hg 18-20 (at the ankle) Antiembolism thigh length with inspection hole and silicon top. Comp: 82% Polyamide, 18% Elastane



Art. H53 - AGT - mm/Hg 18-20 (at the ankle) Antiembolism full length stocking with inspection hole, adjustable top, ambidextrous. Comp: 82% Polyamide, 18% Elastane





CIRCUMFERENCE

NEW		Small	Medium	Large	X-Large	XX-Large
VERSION	<i>l</i> LENGTH	cB 18 - 22 cm cD 29 - 35 cm cG 43 - 54 cm	cB 22 - 26 cm cD 34 - 40 cm cG 48 - 60 cm	cB 26 - 29 cm cD 39 - 45 cm cG 54 - 66 cm	cB 29 - 32 cm cD 44 - 51 cm cG 60 - 72 cm	cB 32 - 35 cm cD 50 - 56 cm cG 67 - 80 cm
Below Knee	Short 33 - 38 cm	H52CS	H52CM	H52CL	H52CXL	H52CXXL
	Normal 38 - 43 cm	H52NS	H52NM	H52NL	H52NXL	H52NXXL
κΒ ∯	Long 43 - 48 cm	H52LS	H52LM	H52LL	H52LXL	H52LXXL
Thigh Length cG	Short 60 - 70 cm	H51CS	H51CM	H51CL	H51CXL	H51CXXL
1 Pair cD	Normal 70 - 80 cm	H51NS	H51NM	H51NL	H51NXL	H51NXXL
¢B ∯	Long 80 - 90 cm	H51LS	H51LM	H51LL	H51LXL	H51LXXL
Full Leg Ambidextrous	Short 60 - 70 cm	H53CS	H53CM	H53CL	H53CXL	H53CXXL
Piece cD /D	Normal 70 - 80 cm	H53NS	H53NM	H53NL	H53NXL	H53NXXL
cB ∯	Long 80 - 90 cm	H53LS	H53LM	H53LL	H53LXL	H53LXXL

